

Protocol
PCR
The Polymerase Chain Reaction

Polymerase Chain Reaction

Amplification of a Short DNA Stretch by Repeated Cycles of In Vitro DNA Polymerization

Science, Vol 239, Issue 4839, 487-491
Copyright © 1988 by American Association for the Advancement of Science

Primer-directed enzymatic amplification of DNA with a thermostable DNA polymerase

RK Saiki, DH Gelfand, S Stoffel, SJ Scharf, R Higuchi, GT Horn, KB Mullis, and HA Erlich
Cetus Corporation, Department of Human Genetics, Emeryville, CA 94608.

A **thermostable DNA polymerase** was used in an **in vitro DNA amplification procedure**, the polymerase chain reaction. The enzyme, isolated from *Thermus aquaticus*, greatly simplifies the procedure and, by enabling the amplification reaction to be performed at higher temperatures, significantly improves the specificity, yield, sensitivity, and length of products that can be amplified. **Single-copy genomic sequences were amplified by a factor of more than 10 million with very high specificity**, and DNA segments up to 2000 base pairs were readily amplified. In addition, the method was used to amplify and detect a target DNA molecule present only once in a sample of 10^5 cells.

- ◆ Mullis, K. B., and F. A. Faloona. "Specific synthesis of DNA in vitro via a [polymerase](#) catalysed chain reaction." *Methods in Enzymology* no. 155 (1987): 335-350.

What is the Polymerase Chain Reaction?

PCR is a method for performing DNA replication in a test tube.

Many, many copies of a small DNA sequence may be amplified from a few starting copies. Once the copies are made, the DNA may more easily be studied. For example, the nucleotide sequence of a particular gene can be determined.

Who uses PCR?:


- PCR is probably the most widely used technique in molecular biology
- This technique is used in biomedical research, criminal forensics, molecular archaeology.

Why is PCR so widely used?:

- It allows you to specifically amplify a single gene or DNA sequence from a very small tissue sample (even a single cell).

Who invented PCR? When was it invented?

- PCR was invented by Kary Mullis (photo right) in 1983.
- Mullis was an employee of the Cetus Corporation.
- He won the Nobel Prize in Chemistry in 1993.
- Since then, he "has become ever more eccentric," according to James D. Watson in *DNA. The Secret of Life*. (*In my opinion, Mr. Watson knows a thing or two about being eccentric.*)




PCR: Polymerase Chain Reaction

How does it work?

- PCR is an *in vitro* process (occurring in a test tube) that imitates the *in vivo* process (occurring in a living cell) of DNA replication.
- There are some similarities and some differences between PCR and DNA replication.

	DNA replication (in vivo)	PCR (in vitro)
What does it accomplish?	Copies the whole genome once.	Makes 2 ³⁰ - 2 ⁴⁰ copies of one gene.
STEPS OF REPLICATION	HOW IS IT ACCOMPLISHED?	HOW IS IT ACCOMPLISHED?
Step 1. Open up double strand of DNA helix.	Enzyme called helicase .	Melt strands apart (heat test tube to near boiling; >90°C)
Step 2. "Prime" the copying process by annealing a small primer to the exposed bases.	Millions of different RNA primers produced by primase - bind throughout the genome.	1 specific pair of DNA primers devised by scientist to match one particular gene.
Step 3. Copy DNA using an enzyme to catalyze the polymerization of individual nucleotide bases.	DNA polymerase + A, C, G, T	DNA polymerase (from hot-spring extremophile). + A, C, G, T.
Is that it?	That's essentially it!	Repeat process 30-40 times, each time doubling the number of copies.

PCR: Polymerase Chain Reaction




Rules of DNA Replication (*in vivo* and *in vitro*)

Base Pairing

- Two simple rules ensure that the newly copied strand has the same sequence as the original strand that it displaced. (1) **A** pairs with **T** via two hydrogen bonds (2) **G** pairs with **C** via three hydrogen bonds.

Directionality of DNA Replication

- Also, remember that DNA replication proceeds from 5' to 3' on each strand, and that the two strands are anti-parallel, so replication is headed in opposite directions on either strand.



Assembling the PCR reaction

For each 20-microliter PCR reaction, in a separate 200-microliter microcentrifuge tube, assemble the following reagents:

Reagents	Stock concentration	Final amount or concentration	Volume per 20-microliter reaction
Reaction buffer	10X	1X	2.0 microliters (μL)
dNTPs	2 mM (10X)	200 μM	2.0 microliters (μL)
thermostable polymerase	5 units / μL	1 unit	0.2 microliters (μL)
forward primer	10 picomoles / μL	20 picomoles	2.0 microliters (μL)
reverse primer	10 picomoles / μL	20 picomoles	2.0 microliters (μL)
template DNA	10-500 nanograms / μL	100-500 ng	1-11.8 microliters (μL)

Troubleshooting PCR

Amplification of the wrong product	Annealing temperature is too low.	Check primer melting temperature. Set annealing temperature to $T_m - 5\text{C}^\circ$
	Mispriming at the correct annealing temperature	Check that primers don't match to multiple sequences in the genome.
	Too much dNTPs	Reduce dNTPs to 25-50 μM
No Product	Annealing temperature is too high.	Check primer melting temperature. Set annealing temperature to $T_m - 5\text{C}^\circ$
	Poor primer design	Verify that primer sequence is correct.
	Missing reagent or not enough of a reagent	Be sure that you are using 20 picomoles of each primer.
	Too few cycles of PCR.	Increase cycle number up to 40.